



Studying the Mechanical Properties of Red Blood Cells with NanoVision

Application Note




Figure 1. A group of red blood cells.




Figure 2. A single red blood cell.




Figure 3. Four indents on a single red blood cell.

Introduction

This application note presents a brief study of the mechanical properties of red blood cells (RBC) using the Agilent NanoVision nanomechanical microscopy option. NanoVision delivers quantitative imaging by coupling a linear electromagnetic actuation-based indentation head with a closed-loop nanopositioning stage on the Agilent G200 Nano Indenter.

Sample Preparation and Imaging

A blood smear was prepared on a microscope slide and was dried for one hour before testing. The microscope slide was then mounted as a regular sample in an Agilent Nano Indenter G200.

Next, the sample is translated in X and Y under a diamond tip. A small load is applied to the tip in order to keep it in contact with the surface. The height variations recorded by the tip create the topographic image of the sample's surface. The piezo-actuated sample stage used in the NanoVision option provides high accuracy in lateral positioning (2 nm) and imaging of the samples. A cube-corner diamond tip is utilized to image the surface profile of the sample.

To provide a high degree of accuracy as well as excellent control for the application of the load during imaging, the Agilent Dynamic Contact Module (DCM) indentation head option is also used with the G200.

The image in Figure 1 shows the well-known biconcave disk shape of the red blood cells. A single cell (Figure 2) was selected for the first study of mechanical properties.

Indentation

Indentation experiments were performed using Continuous Stiffness Measurement (CSM). With the option elastic modulus and hardness are measured as a continuous function of penetration. (Without this option measurement of elastic modulus and hardness can only be achieved at the maximum penetration depth).

Two series of indentations were performed on the red blood cell sample. The first group of four indents was performed on a single red blood cell (Figure 3). The four indents were made at a constant strain rate of 0.2/sec and with a maximum penetration depth of 1000 nm. The hardness and modulus measurements

(Figures 4 and 5, respectively) present very good repeatability. The substrate effect is detected around 100nm on the modulus.

The second group of indents was made on several cells. One indent was made at the center of each cell. The eight indents were made at a constant strain rate of 0.2/sec and with a maximum penetration depth of 500nm. The results for the measurements of hardness and modulus are consistent from cell to cell (Figures 6 and 7, respectively). The indents were imaged via the same technique used previously.

Conclusions

The use of NanoVision in conjunction with a Nano Indenter G200 allows both indentation and imaging of biological materials, such as red blood cells. Measurement of mechanical properties of such materials using this instrumentation yields highly accurate and repeatable data.

Technology and Applications

The Nano Indenter G200 is powered by electromagnetic actuation to achieve unparalleled dynamic range in force and displacement. The instrument's unique design avoids lateral displacement artifacts, while software compensates fully for any drift in force. The DCMII option offers impressive performance and advantages including 3x higher loading capability (30 mN max load), easy tip exchange for quick removal and installation of application-specific tips, and a full 70 μm range of indenter travel. Using the G200, researchers can measure Young's modulus and hardness in compliance with ISO 14577. Deformation can be measured over six orders of magnitude (from nanometers to millimeters).

Applications include semiconductor, thin films, and MEMs (wafer applications); hard coatings and DLC films; composite materials, fibers, and polymers; metals and ceramics; and biomaterials and biology.




Figure 4. Hardness for the four indents in a single RBC.




Figure 5. Modulus for the four indents in a single RBC.




Figure 6. Average of Hardness for eight indents on different RBC.




Figure 7. Average of Modulus for eight indents on different RBC.

Nano Mechanical Systems from Agilent Technologies

Agilent Technologies, the premier measurement company, offers high-precision, modular nano-measurement solutions for research, industry, and education. Exceptional worldwide support is provided by experienced application scientists and technical service personnel. Agilent's leading-edge R&D laboratories ensure the continued, timely introduction and optimization of innovative, easy-to-use nanomechanical system technologies.

www.agilent.com/find/nanoindenter

Americas

Canada	(877) 894 4414
Latin America	305 269 7500
United States	(800) 829 4444

Asia Pacific

Australia	1 800 629 485
China	800 810 0189
Hong Kong	800 938 693
India	1 800 112 929
Japan	0120 (421) 345
Korea	080 769 0800
Malaysia	1 800 888 848
Singapore	1 800 375 8100
Taiwan	0800 047 866
Thailand	1 800 226 008

Europe & Middle East

Austria	43 (0) 1 360 277 1571
Belgium	32 (0) 2 404 93 40
Denmark	45 70 13 15 15
Finland	358 (0) 10 855 2100
France	0825 010 700*
	*0.125 €/minute
Germany	49 (0) 7031 464 6333
Ireland	1890 924 204
Israel	972-3-9288-504/544
Italy	39 02 92 60 8484
Netherlands	31 (0) 20 547 2111
Spain	34 (91) 631 3300
Sweden	0200-88 22 55
Switzerland	0800 80 53 53
United Kingdom	44 (0) 118 9276201

Other European Countries:
www.agilent.com/find/contactus

Product specifications and descriptions in this document subject to change without notice.

© Agilent Technologies, Inc. 2009
Printed in USA, November 12, 2009
5990-5020EN